




# YOUR WAY TO FITTING EXCELLENCE




AURICAL® Aud



otometrics

# THE AUDIOMETER

655g of concentrated knowledge



**Unique design:** Compact and designed to support your needs.



**Connectivity:** The vast connectivity possibilities make it the natural central point of the fitting solution.

**Under desk mounting:** Mounting the audiometer with the built-in HI-PRO 2 under the desk helps tidy up cabling and provides convenient access to the HI-PRO plugs.



See the videos about the complete fitting solution at [www.otometrics.com/aurical](http://www.otometrics.com/aurical)



AURICAL Aud is the world's only true fitting audiometer. It integrates an audiometer, counseling tools, power amplifier, USB hub and optional HI-PRO 2 with High Speed Mode in a modern, ultra-compact design. It's a complete, portable audiometer that offers outstanding flexibility. It is designed to fit into all kinds of environments, inside and outside your clinic.

The audiometer is easy to use too. The exceptionally intuitive OTOSuite® interface has an on-screen control panel and easy-to-access User Test functions that can be customized to suit individual needs. It ensures fast operation and consistently accurate results.

High quality sound is important for your fitting results. The unique RoomTune™ functionality turns AURICAL Aud into a compensated sound system allowing you to present more accurate and reliable sound from your fitting modules.

## THE FACTS

- AC, BC and SF audiometry (built-in power amplifier)
- Tone audiometry testing up to 12.5 kHz, including FRESH noise™
- Loudness Scaling (optional)
- Integrated speech material
- SAL, TEN, QuickSIN tests (some optional)
- Masking assistant
- Sound adjustment with RoomTune™
- Built-in USB powered
- Optional HI-PRO 2 w/ High Speed Mode
- OTOSuite with User Tests
- Intuitive on-screen control panel
- One-click data logging
- Selection of report templates



**Quick Select**

AC: [R] [B] [L] [X]

Mask: [A] [U]

BC: [K] [A] [G] [J]

Mask: [E] [N]

SF: [S] [S] [S]

Mask: [X] [X]

Condition: Unaided [Threshold] [MCL] [UCL]

**Test Options**

**dB Step**: 1 [5] [10]

**AC Transducer**: [Insert] [Phone]

**Control**: [Stim. On] [Stim. Lock] [Tracking] [Pulsed]

**Stimulus**: [Tone] [Warble] [Stimulus 2]

**Monitor / Level**

**MONITOR**: [Stim] [Mask] [Talk back]

**LEVEL**: [Mic] [Source A] [Source B]

**Quick Select**

**Test Ear**: [R] [B] [L]

**Measurement**: [AC] [BC] [SF] [SDT] [SRT] [MCL] [UCL] [WRS/SRS] [-SNR] [Off]

[Unaided] [Ipsi] [Masking On]

**Test Options**

**dB Step**: 1 [5] [10]

**AC Transducer**: [Insert] [Phone]

**Control**: [Stim. On] [Stim. Lock] [Tracking]

**Input**: [Int. File] [Source A]

**Monitor / Level**

**MONITOR**: [Stim] [Mask] [Talk back]

**LEVEL**: [Mic] [Source A] [Source B]

**Test Selector**

**User Test:**

- Adult - Extensive Aud
- Adult - Initial Aud
- HI Benefit Aided Speech SF
- HI Benefit Aided Speech SF in Noise
- Pediatric - Extensive Aud
- Pediatric - Initial Aud

**Special Test:**

1. Full Data Overview
2. ABLB
3. SISI
4. Weber (+ Rinne Data Entry)
5. Stenger
6. Tone Decay
- Default USA

[Cancel] [Help]

Various view options allow you to choose what is best for you and your patient. Go from simple to advanced view in 2 mouse clicks.

The on-screen control panels in AURICAL Aud and OTOsuite are intuitive and fast.

# THE AURICAL FAMILY

Fully integrated with OTOsuite



The OTOsuite® software integrates all the elements of AURICAL®. It's intuitive, easy to learn and easy to use. It supports the entire client journey, from diagnostics and counseling, to fitting and verification. Simple to use user test functionality ensures configurable, consistent and smooth workflows that can be customized to suit individual needs. And you work in a single software environment across modules. The result is a better fitting experience – for you and for your client.

## About Otometrics

Otometrics is the audiology industry leader providing instrumentation and software solutions to hearing and balance care professionals worldwide. For more than 50 years, Otometrics has been helping hearing and balance care professionals succeed in improving the quality of life for their clients and patients around the world by delivering expert knowledge, reliable solutions and trusted partnership.

Otometrics develops, manufactures and markets computer-based audiological, otoneurologic and vestibular instrumentation in more than 85 countries under the MADSEN®, AURICAL®, HORTMANN® and ICS® brand names. Otometrics is part of the GN Group, one of the largest companies in Denmark. For more information, please visit [www.otometrics.com](http://www.otometrics.com).