



Amtronix
Diagnostics

Lebone Medical Supplies

RT-274 2023 CATALOGUE

Audiological equipment,
consumables and bone
conduction hearing aids

www.amtronixdiagnostics.co.za



Screening audiometers

Spec sheet:



RT274-01-022 OPTION 1

R48 914,10

INVENTIS TRIANGLE AUDIOMETER

- 125-8000 Hz one-channel portable audiometer
- Battery / plug-in operated
- Unmasked puretone air and bone conduction testing
- Transducers: earphones, bone conductor, response button
- Carrycase included
- AC: -10 to 100 dB HL
- BC: -10 to 60 dB HL

Spec sheet:



RT274-01-022 OPTION 2

R56 315,50

MAICO MA28 AUDIOMETER

- 125-8000 Hz One-channel portable audiometer
- Plug-in operated
- Unmasked puretone air and bone conduction testing
- Transducers: earphones, bone conductor, response button
- Carrycase included
- AC: -10 dB HL to 110 dB HL
- BC: -10 dB HL to 80 dB HL

Spec sheet:



RT274-01-023 OPTION 1

R87 204,50

PATH SENTI FLEX HF AUDIOMETER

- 125-16 000Hz One-channel portable audiometer
- Battery or plug-in operated
- Unmasked puretone air and bone conduction testing
- Transducers: earphones, bone conductor, response button
- Carrycase included

(Additional capabilities with add-on software and transducers: boothless audiometry, recorded speech & MAGIC (play audiometry))



Diagnostic audiometers

INVENTIS HARP PLUS



RT274-01-050

R130 628,50

INVENTIS HARP PLUS DIAGNOSTICS AUDIOMETER WITH LAPTOP

RT274-01-056

R118 553,50

INVENTIS HARP PLUS DIAGNOSTICS AUDIOMETER ONLY

- Two-channel audiometer
- Masked puretone and warbletone air and bone conduction and speech testing
- Accessories: Earphones, insert phones, bone conductor, response button, talk-forward, built-in or gooseneck microphone
- Pre-recorded speech material in South African accents and languages
- USB ports to connect to an external CD player for recorded tests
- Special tests: Auto threshold ,SISI, ABLB, Stenger, Bekesy, Weber, tone decay, DLI
- 125-8 000Hz from -10-120dB in increments of 5db.

SPEC SHEET:



INVENTIS PIANO HF



RT274-01-051

R188 186,00

INVENTIS PIANO CLINICAL HF AUDIOMETER WITH LAPTOP

RT274-01-057

R176 111,00

INVENTIS PIANO CLINICAL HF AUDIOMETER ONLY

- Two-channel audiometer
- Masked puretone and warbletone air and bone conduction and speech testing
- Accessories: DD45 Earphones, HDA-200 earphones, insert phones, bone conductor, response button, talk-forward, built-in or gooseneck microphone
- Pre-recorded speech material in South African accents and languages
- USB ports to connect to an external CD player for recorded tests,
- Special tests: Auto threshold ,SISI, ABLB, Stenger, Bekesy, Weber, tone decay, DLI, TEN test.
- 125-20 000Hz from -10-120dB in increments of 5db.

SPEC SHEET:





Screening tympanometers

MAICO EASYTYMP: SCREENING TYMP & IPSI REFLEXES



RT274-01-024 OPTION ONE

R78 619,29

- 226Hz probe tone
- Pressure: -400-+200daPa
- Ipsilateral reflexes: 500,1000,2000,4000Hz
- Automated screening
- Full colour high resolution display
- Stores more than 1000 measurements
- Battery operated
- Included: Thermal printer with 2 rolls of paper, carry case, eartip box, calibration test cavity

Optional add-ons: PC software, 1000Hz probetone, BBN reflex, cradle kit

SPEC SHEET:



INVENTIS TIMPANI: SCREENING TYMP & IPSI REFLEXES



RT274-01-024 OPTION TWO

R78 081,55

- 226Hz probe tone
- Pressure: -400- +200daPa
- Ipsilateral reflexes: 500,1000,2000,4000Hz
- Automated screening
- Full colour high resolution display
- Stores more than 1000 measurements
- Battery operated
- Included: Thermal printer with 2 rolls of paper, carry case, eartip box, calibration test cavity

Optional add-ons: 1000Hz probetone, screening puretone audiometry, PC software

SPEC SHEET:





Diagnostic Tympanometers

INVENTIS CLARINET

MAICO MI34 TOUCHTYMP



RT274-01-052 OPTION ONE

R151 926,50

INVENTIS CLARINET BASIC INCLUDING A LAPTOP AND SOFTWARE

RT274-01-058 OPTION ONE

R136 401,50

INVENTIS CLARINET BASIC ONLY

RT274-01-059 OPTION ONE

R148 292,50

INVENTIS CLARINET PLUS ONLY

RT274-01-053 OPTION ONE

R163 817,50

INVENTIS CLARINET PLUS INCLUDING A LAPTOP AND SOFTWARE

RT274-01-052 OPTION TWO

R150 712,10

MAICO MI34 BASIC INCLUDING A LAPTOP AND SOFTWARE

RT274-01-058 OPTION TWO

R133 462,10

MAICO MI34 BASIC ONLY

RT274-01-059 OPTION TWO

R143 835,10

MAICO MI34 HF ONLY

RT274-01-053 OPTION TWO

R161 085,10

MAICO MI34 HF INCLUDING A LAPTOP AND SOFTWARE

INVENTIS CLARINET BASIC

- Diagnostic tympanometry (226Hz probe tone)
- Pressure: -600- +400
- ipsi- and contra-lateral reflexes (500, 1000, 2000, 4000Hz, BBN, LPN, HPN)
- Reflex decay
- Eustachian tube function (intact and perforated tympanic membrane),
- USB PC interface.
- Included: Eartip box and calibration cavity

INVENTIS CLARINET PLUS

Includes above features plus

- 678Hz, 800Hz and 1000Hz probe tone
- Multicomponent tympanometry



SPEC SHEET:

MAICO MI34 BASIC

- Diagnostic tympanometry (226, 678,800Hz Probe tone)
- Pressure: -600- +400
- ipsi- and contra-lateral reflexes (500, 1000, 2000, 4000Hz, BBN, LPN, HPN)
- Reflex decay
- Eustachian tube function (intact and perforated tympanic membrane),
- USB PC interface.
- Included: Eartip box and calibration cavity

MAICO MI34 BASIC HF

Includes above features plus

- 1000Hz probe tone



SPEC SHEET:



OAE and ABR

Spec sheet:



RT274-01-054

(INCLUDES OAE+ LAPTOP)

R139 092,50

RT274-01-060

(INCLUDES OAE ONLY)

R121 842.50

PATH DIAGNOSTIC OAE

- Diagnostic high frequency DPOAE and TEOAE base unit:
- Six configurable protocols testing from
- Stores >500 tests
- Displays normative data
- includes: OAE probe, probe calibration test cavity, Mira database software, micro USB charging / transfer cable, carry case, eartip kit and replacement probe tubes I and calibration certificate.

(Additional capabilities with add-on software and transducers: audiometry, recorded speech, MAGIC (play audiometry))

Spec sheet:



RT274-01-025

OPTION ONE

R127 339,50

PATH ADVANCED OAE + AABR

- Screening DPOAE + AABR
- Built in normative values
- Includes: USB, software, OAE probe, ABR cable with 3 electrode leads, charger, carrycase, probetip and eartip box
- Stores >500 tests

Additional capabilities with add-on software and transducers: audiometry, MAGIC (play audiometry), recorded speech, OAEs, Immittance, AABR, ABR, ASSR, ECOCHG, VEMP)

Spec sheet:



RT274-01-025

OPTION TWO

R199 582,50

MAICO EASYSCREEN OAE + AABR (Beraphone)

- Screening DPOAE + AABR
- Built in normative values
- Includes: USB, software, OAE probe, ABR cable and beraphone, charger, carrycase, probetip and eartip box
- Stores >500 tests



Path Medical Equipment



**PATH SENTIERO ADVANCED
CHASSIS ONLY
(NO MODULES)**

RT274-01-036

R72 392,50

PATH SENTIERO ADVANCED CHASSIS PLUS ONE MODULE

Path Sentiero Advanced with screening AC+BC audiometer module

RT274-01-037 R111 193,50

Path Sentiero Advanced with HF screening AC+BC audiometer module

RT274-01-038 R137 758,50

Path Sentiero Advanced with screening 226Hz tymps+ reflexes module

RT274-01-039 R111 550,00

Path Sentiero Advanced with 226, 678+100Hz tymps+ reflexes module

RT274-01-040 R125 465,00

Path Sentiero Advanced with screening DPOAE & TEOAE module

RT274-01-041 R126 026,20





Path Sentiero Advanced with high frequency DPOAE module

RT274-01-042 R127 443,00

Path Sentiero Advanced with diagnostic ASSR module

RT274-01-063 R151 627,50

ADD-ON MODULE FOR PATH MEDICAL DEVICES

	SENTIERO Advanced	SENTIERO Diagnostic	SENTIERO Desktop	SENTI Desktop	
ADD-ON MODULE					
AUDIO screening AC+BC	●	●	●	●	RT274-01-043 R49 093,50
HF AUDIO screening AC+BC	●	●	●	●	RT274-01-044 R75 658,50
DPOAE+TEOAE screening	●	●	●		RT274-01-047 R55 761,20
DPOAE Diagnostic	●	●	●		RT274-01-048 R59 834,50
Tympanometry screening 226Hz tymps+ reflexes	add on Tymp		●		RT274-01-045 R53 475,00
Tympanometry 226, 678+100Hz tymps+ reflexes	add on Tymp		●		RT274-01-046 R67 459,00
AABR	●				RT274-01-049 R74 669,50

Note: above codes include the test's respective transducers and probe tips



Real Ear Measures. RECD. Hearing Instrument test box.

INVENTIS TRUMPET AND DRUM REM AND HIT SYSTEM

RT274-01-064

R203 653,50



Trumpet
Spec sheet:



Drum
Spec sheet:



REUR / REUG, REAR / REAG, REIG, REOR / REOG, MPO, Advanced, Live speech, RECD,
Paired insert phones
in-situ headset
probe kit (tubes, O-rings)
Monitor headset
Software USB
Medical grade power supply
USB connection cable

Additional add-on capabilities: Puretone audiometry, speech audiometry.



Probe tip and ear tips kits

PATH MEDICAL Paediatric Ear Tip Kit (Tymps+ OAEs)




Path Medical Newborn and Baby Ear Tips, Probe Tips and Test Cavity

- Path Medical Ear Tip Kit**
- 2 Paediatric probe tips
 - 1 probe test cavity
 - 120 ear tips

RT274-01-004/5
R3 522,45

PATH MEDICAL Paediatric, Young children and Adult Ear Tip Kit (Tymps+ OAEs)



NEWBORN

CHILD/ADULT

- Path Medical Ear Tip Kit**
- 4 Paediatric probe tips
 - 4 Adult probe tips
 - 2 probe test cavities
 - 80 ear tips

RT274-01-006
R3 522,45

Probe tip and ear tips kits

NB: when ordering with the code : RT274-01-007/8, specify whether you would like the Inventis or the Maico eartip box



inventis
Audiology Equipment

 WHITE Ø 6 mm IN/C1124	 BLUE Ø 14 mm IN/C1129	 ORANGE Ø 15 mm IN/C1130	 LIGHT BLUE Ø 18 mm IN/C1131
 inventis Audiology Equipment			
 PINK Ø 7 mm IN/C1125	 GREEN Ø 10 mm IN/C1127	 RED Ø 12 mm IN/C1128	 YELLOW Ø 21 mm IN/C1132

(Note: The image also shows a photograph of the eartips in their respective compartments within the kit box.)

Inventis Eartip kit

- Suitable for the Inventis tympanometers
- Contains eartips in 8 variations
- Suitable eartips for adults and children
- 1 set contains 160 pieces

RT274-01-007/8

R3 775,45



MAICO

(Note: The image shows an open Maico Eartip kit box containing various colored eartips and cleaning tools.)

Maico Eartip kit

- Suitable for the Maico EasyTymp, Maico touchTymp
- Contains eartips in 9 variations
- Suitable eartips for adults and children
- Including equipment to remove the tips from the tympanometer
- Including ear loop to clean the eartips if needed
- 1 set contains 70 pieces

• NB: specify

RT274-01-007/8

R3 775,45



Audiology tools and consumables




RT274-02-009

HEINE SET OF 4 SPECULAE.
2.4MM, 3.0MM, 4.0MM, 5.0 M
R994,75



RT274-02-001

OTOLIGHT
R517,50



RT274-02-004

REUSABLE OTOLIGHT TIPS X2
R256,45



RT274-07-048

PURESIL IMPRESSION MATERIAL 1:1
R1 159,20



RT274-07-049

PURESIL IMPRESSION MATERIAL
R1 159,20



RT274-07-050

PURESIL IMPRESSION CATALYST
R454,25



RT274-07-043

PURETONE OTOSTOPS
SMALL
R359,95



RT274-07-044

PURETONE OTOSTOPS
MEDIUM
R359,95



RT274-07-045

PURETONE OTOSTOPS
LARGE
R359,95



Audiology tools and consumables



RT274-01-013

EEG ABRASIVE SKIN PREPPING GEL
R465,75



RT274-01-012

EEG TEN 20 ELECTRODE PASTE
OR SONAVELLE ELECTRODE GEL
(SPECIFY WHICH ONE YOU WANT)
R2 202,25



RT274-02-012

AUDIOMETER EARPHONE COVERS, ,
ACOUSTICALLY TRANSPARENT.
PACK OF 50.
R466,90




RT274-01-020

1X REUSABLE GOLD CUP
ELECTRODE LEAD
ADULT-10MM
R448,50



RT274-01-021

1X REUSABLE GOLD CUP
ELECTRODE LEAD
PEADIATRIC-6MM
R448,50



RT274-01-011

INSERT FOAM EARTIPS
MEDIUM-13.7MM (50 PER PACK) **OR**
OR SMALL-9.4MM (50 PER PACK)
SPECIFY SIZE WANTED
R859,05



RT274-03-013

HARTMAN CROCODILE FORCEPS
R966,00



RT274-03-007

SYRINGE, EAR IRRIGATION, PLASTIC
R189,75

Note: product supplied may not look exactly like the picture (EG: colour may differ depending on supplier stock)



Hearing Aid Care



RT274-07-021

PURETONE THIN-WALLED TUBING
R399,05



RT274-07-022

PURETONE MEDIUM-WALLED TUBING
R399,05



RT274-07-023

PURETONE THICK-WALLED TUBING
R399,05



RT274-02-014

PURETONE BLOWER
R281,75



RT274-02-008

STETOCLIP
R937,25



RT274-07-036

PAEDIATRIC HEARING AID CLIP
R327,75



RT274-07-013

BATTERY TESTER
R238,05



RT274-07-034

HYGIENE KIT
CONTAINS 6 PURECLEAN CLEANING WIPES, A PURECLEAN CLEANING SPRAY, A PURECLEAN DRYING CAPSULE AND A PURECLEAN DRYING CUP. INCLUDES ZIP BAG, BROCHURE AND PURECLEAN HEARING AIDS SLEEVE.
R1519,15



RT274-07-040

DRY AID, ELECTRONIC, FOR HEARING AIDS, UNIVERSAL, KIT
R2 340,25

Note: product supplied may not look exactly like the picture (EG: colour may differ depending on supplier stock)



Bone conduction hearing aids

Spec sheet:



PURETONE SOFTBAND MONAURAL BONE CONDUCTION HEARING AID

For mild-moderate conductive loss
Softband colour options: Black, blue, red or white
(Output 30-40)
6 Batteries (Size 13)
R15,410.00

**RT274-04-004
OFFER 2**

Spec sheet:



PURETONE DIGITAL HEADBAND MONAURAL BONE CONDUCTION HEARING AID

6 BATTERIES (SIZE 13)
R16 675

**RT274-04-005
OFFER 2**

**Audiology
equipment and
consumables**



Lebone Medical Supplies

Trading as

Amtronix Diagnostics

Partnering with you.



Front desk

Name. Thompho Lumadi
 Phone. 012 054 5100/ 073 823 2380
 Email. hear4u@amtronix.co.za

Audiology: Support, equipment advice, training

Name. Kayla Mc Farlane/ Fehr Tshiteya
 Phone. 081 046 1626/ 078 498 3028
 Email. audiochat@amtronix.co.za / listen@amtronix.co.za

Sales, quotes & orders

Name. Merle Pemberton
 Phone. 078 350 5483
 Email. diagnostic@amtronix.co.za

Arrangement of calibrations, call-outs & repairs

Name: Sorette Olivier
 Phone. 012 054 5100/ 073 823 2380
 Email. allears@amtronix.co.za

Technical support

Name. Tiaan Schoeman
 Phone. 060 558 6199
 Email. 911@amtronix.co.za



Website



Most up-to date
Catalogue of audiology
equipment on tender