



Amtronix
Diagnostics

Lebone Medical Supplies

RT-274 2023 CATALOGUE

Audiological equipment,
consumables and bone
conduction hearing aids

www.amtronixdiagnostics.co.za



Screening audiometers

Spec sheet:



INVENTIS TRIANGLE AUDIOMETER

- 125-8000 Hz one-channel portable audiometer
- Battery / plug-in operated
- Unmasked puretone air and bone conduction testing
- Transducers: earphones, bone conductor, response button
- Carrycase included
- AC: -10 to 100 dB HL
- BC: -10 to 60 dB HL

RT274-01-022 OPTION 1

R48 914,10

Spec sheet:



MAICO MA28 AUDIOMETER

- 125-8000 Hz One-channel portable audiometer
- Plug-in operated
- Unmasked puretone air and bone conduction testing
- Transducers: earphones, bone conductor, response button
- Carrycase included
- AC: -10 dB HL to 110 dB HL
- BC: -10 dB HL to 80 dB HL

RT274-01-022 OPTION 2

R56 315,50

Spec sheet:



PATH SENTI FLEX HF AUDIOMETER

- 125-16 000Hz One-channel portable audiometer
- Battery or plug-in operated
- Unmasked puretone air and bone conduction testing
- Transducers: earphones, bone conductor, response button
- Carrycase included

(Additional capabilities with add-on software and transducers: boothless audiometry, recorded speech & MAGIC (play audiometry))

RT274-01-023 OPTION 1

R87 204,50



Diagnostic audiometers

INVENTIS HARP PLUS



RT274-01-050

R130 628,50

INVENTIS HARP PLUS DIAGNOSTICS AUDIOMETER WITH LAPTOP

RT274-01-056

R118 553,50

INVENTIS HARP PLUS DIAGNOSTICS AUDIOMETER ONLY

- Two-channel audiometer
- Masked puretone and warbletone air and bone conduction and speech testing
- Accessories: Earphones, insert phones, bone conductor, response button, talk-forward, built-in or gooseneck microphone
- Pre-recorded speech material in South African accents and languages
- USB ports to connect to an external CD player for recorded tests
- Special tests: Auto threshold ,SISI, ABLB, Stenger, Bekesy, Weber, tone decay, DLI
- 125-8 000Hz from -10-120dB in increments of 5db.

SPEC SHEET:



INVENTIS PIANO HF



RT274-01-051

R188 186,00

INVENTIS PIANO CLINICAL HF AUDIOMETER WITH LAPTOP

RT274-01-057

R176 111,00

INVENTIS PIANO CLINICAL HF AUDIOMETER ONLY

- Two-channel audiometer
- Masked puretone and warbletone air and bone conduction and speech testing
- Accessories: DD45 Earphones, HDA-200 earphones, insert phones, bone conductor, response button, talk-forward, built-in or gooseneck microphone
- Pre-recorded speech material in South African accents and languages
- USB ports to connect to an external CD player for recorded tests,
- Special tests: Auto threshold ,SISI, ABLB, Stenger, Bekesy, Weber, tone decay, DLI, TEN test.
- 125-20 000Hz from -10-120dB in increments of 5db.

SPEC SHEET:





Screening tympanometers

MAICO EASYTYMP: SCREENING TYMP & IPSI REFLEXES



RT274-01-024 OPTION ONE

R78 619,29

- 226Hz probe tone
- Pressure: -400-+200daPa
- Ipsilateral reflexes: 500,1000,2000,4000Hz
- Automated screening
- Full colour high resolution display
- Stores more than 1000 measurements
- Battery operated
- Included: Thermal printer with 2 rolls of paper, carry case, eartip box, calibration test cavity

Optional add-ons: PC software, 1000Hz probetone, BBN reflex, cradle kit

SPEC SHEET:



INVENTIS TIMPANI: SCREENING TYMP & IPSI REFLEXES



RT274-01-024 OPTION TWO

R78 081,55

- 226Hz probe tone
- Pressure: -400- +200daPa
- Ipsilateral reflexes: 500,1000,2000,4000Hz
- Automated screening
- Full colour high resolution display
- Stores more than 1000 measurements
- Battery operated
- Included: Thermal printer with 2 rolls of paper, carry case, eartip box, calibration test cavity

Optional add-ons: 1000Hz probetone, screening puretone audiometry, PC software

SPEC SHEET:





Diagnostic Tympanometers

INVENTIS CLARINET



RT274-01-052 OPTION ONE

R151 926,50

INVENTIS CLARINET BASIC INCLUDING A LAPTOP AND SOFTWARE

RT274-01-058 OPTION ONE

R136 401,50

INVENTIS CLARINET BASIC ONLY

RT274-01-059 OPTION ONE

R148 292,50

INVENTIS CLARINET PLUS ONLY

RT274-01-053 OPTION ONE

R163 817,50

INVENTIS CLARINET PLUS INCLUDING A LAPTOP AND SOFTWARE

INVENTIS CLARINET BASIC

- Diagnostic tympanometry (226Hz probe tone)
- Pressure: -600- +400
- ipsi- and contra-lateral reflexes (500, 1000, 2000, 4000Hz, BBN, LPN, HPN)
- Reflex decay
- Eustachian tube function (intact and perforated tympanic membrane),
- USB PC interface.
- Included: Eartip box and calibration cavity

INVENTIS CLARINET PLUS

Includes above features plus

- 678Hz, 800Hz and 1000Hz probe tone
- Multicomponent tympanometry



SPEC SHEET:

MAICO MI34 TOUCHTYMP



RT274-01-052 OPTION TWO

R150 712,10

MAICO MI34 BASIC INCLUDING A LAPTOP AND SOFTWARE

RT274-01-058 OPTION TWO

R133 462,10

MAICO MI34 BASIC ONLY

RT274-01-059 OPTION TWO

R143 835,10

MAICO MI34 HF ONLY

RT274-01-053 OPTION TWO

R161 085,10

MAICO MI34 HF INCLUDING A LAPTOP AND SOFTWARE

MAICO MI34 BASIC

- Diagnostic tympanometry (226, 678,800Hz Probe tone)
- Pressure: -600- +400
- ipsi- and contra-lateral reflexes (500, 1000, 2000, 4000Hz, BBN, LPN, HPN)
- Reflex decay
- Eustachian tube function (intact and perforated tympanic membrane),
- USB PC interface.
- Included: Eartip box and calibration cavity

MAICO MI34 BASIC HF

Includes above features plus

- 1000Hz probe tone



SPEC SHEET:



OAE and ABR

Spec sheet:



RT274-01-054

(INCLUDES OAE+ LAPTOP)

R139 092,50

RT274-01-060

(INCLUDES OAE ONLY)

R121 842.50

PATH DIAGNOSTIC OAE

- Diagnostic high frequency DPOAE and TEOAE base unit:
- Six configurable protocols testing from
- Stores >500 tests
- Displays normative data
- includes: OAE probe, probe calibration test cavity, Mira database software, micro USB charging / transfer cable, carry case, eartip kit and replacement probe tubes I and calibration certificate.

(Additional capabilities with add-on software and transducers: audiometry, recorded speech, MAGIC (play audiometry))

Spec sheet:



RT274-01-025

OPTION ONE

R127 339,50

PATH ADVANCED OAE + AABR

- Screening DPOAE + AABR
- Built in normative values
- Includes: USB, software, OAE probe, ABR cable with 3 electrode leads, charger, carrycase, probetip and eartip box
- Stores >500 tests

Additional capabilities with add-on software and transducers: audiometry, MAGIC (play audiometry), recorded speech, OAEs, Immittance, AABR, ABR, ASSR, ECOCHG, VEMP)

Spec sheet:



RT274-01-025

OPTION TWO

R199 582,50

MAICO EASYSscreen OAE + AABR (Beraphone)

- Screening DPOAE + AABR
- Built in normative values
- Includes: USB, software, OAE probe, ABR cable and beraphone, charger, carrycase, probetip and eartip box
- Stores >500 tests



Path Medical Equipment



**PATH SENTIERO ADVANCED
CHASSIS ONLY
(NO MODULES)**

RT274-01-036

R72 392,50

PATH SENTIERO ADVANCED CHASSIS PLUS ONE MODULE

Path Sentiero Advanced with screening AC+BC audiometer module

RT274-01-037 R111 193,50

Path Sentiero Advanced with HF screening AC+BC audiometer module

RT274-01-038 R137 758,50

Path Sentiero Advanced with screening 226Hz tymps+ reflexes module

RT274-01-039 R111 550,00

Path Sentiero Advanced with 226, 678+1000Hz tymps+ reflexes module

RT274-01-040 R125 465,00

Path Sentiero Advanced with screening DPOAE & TEOAE module

RT274-01-041 R126 026,20

Path Sentiero Advanced with high frequency DPOAE module

RT274-01-042 R127 443,00

Path Sentiero Advanced with diagnostic ASSR module

RT274-01-063 R151 627,50

ADD-ON MODULE FOR PATH MEDICAL DEVICES

COMPARISON CHART

ADD-ON MODULE	SENTIERO Advance	SENTIERO Diagnostic	SENTIERO Desktop	SENTI Desktop	
AUDIO screening AC + BC	✓	✓	✓	✓	RT274-01-043 R49 093,50
HF AUDIO SCREENING AC + BC	✓	✓	✓	✓	RT274-01-044 R75 658,50
DPOAE + TEOAE screening	✓	✓	✓	✗	RT274-01-047 R55 761,20
DPOAE Diagnostic	✓	✓	✓	✗	RT274-01-048 R59 834,50
Tympanometry screening 226Hz tymps + reflexes	add on Tymp	✗	✓	✗	RT274-01-045 R53 475,00
Tympanometry screening 226, 678 + 1000Hz tym + reflexes	add on Tymp	✗	✓	✗	RT274-01-046 R67 459,00
AABR	✓	✗	✗	✗	RT274-01-049 R74 669,50

Note: above codes include the test's respective transducers and probe tips



Real Ear Measures. RECD. Hearing Instrument test box.

INVENTIS TRUMPET AND DRUM REM AND HIT SYSTEM



RT274-01-064

R203 653,50

Trumpet
Spec sheet:



Drum
Spec sheet:



REUR / REUG, REAR / REAG, REIG, REOR / REOG, MPO, Advanced, Live speech, RECD,
Paired insert phones
in-situ headset
probe kit (tubes, O-rings)
Monitor headset
Software USB
Medical grade power supply
USB connection cable

Additional add-on capabilities: Puretone audiometry, speech audiometry.



Probe tip and ear tips kits

PATH MEDICAL Paediatric Ear Tip Kit (Tymps+ OAEs)

5x Path Medical white probetip (newborn)

1x Path Medical Test Cavity for Probe Tip (Newborn)- Blue

20x Path Medical Ear Tips 3.7mm Transparent

20x Path Medical Ear Tips 4mm Red

20x Path Medical Ear Tips 4.5mm Blue

20x Path Medical Ear Tips 5mm (Dark Blue)

20x Path Medical Ear Tips 4.6-7mm (Blue)

20x Path Medical Ear Tip Soft 8mm- broad range

Path Medical Newborn and Baby Ear Tips, Probe Tips and Test Cavity

Neonates - Path Medical Ear Tip Kit

- 2 Paediatric probe tips
- 1 probe test cavity
- 120 ear tips

RT274-01-004

R3 522,45

Paediatrics - Path Medical Ear Tip Kit

- 2 Paediatric probe tips
- 1 probe test cavity
- 120 ear tips

RT274-01-005

R3 522,45

PATH MEDICAL Paediatric, Young children and Adult Ear Tip Kit (Tymps+ OAEs)

White probe tip (newborn)

Test Cavity for Probe Tip (Newborn)

Clear probe tip (Adult)

Test Cavity for Probe Tip (Adult)

Ear Tips 4mm Red

Ear Tips 4.5mm Blue

Ear Tips 4.6-7mm (Blue)

Ear Tips Lamella 8-11mm (Transparent)

Ear Tips soft 6-8mm

Ear Tips 5mm (Dark Blue)

Ear Tips Soft 14mm

Ear Tips soft (10mm)

NEWBORN

Lamella Ear Tips 13-16mm (Blue)

Ear Tip soft 18mm

CHILD/ADULT

Adults - Path Medical Ear Tip Kit

- Paediatric probe tips
- Adult probe tips
- Probe test cavities
- 80 ear tips

RT274-01-006

R3 522,45

Note: product supplied may not look exactly like the picture (EG: colour may differ depending on supplier stock)



Probe tip and ear tips kits

NB: when ordering with the code : RT274-01-007/8, specify whether you would like the Inventis or the Maico eartip box



Inventis Eartip kit

- Suitable for the Inventis tympanometers
- Contains eartips in 8 variations
- Suitable earpieces for adults and children
- 1 set contains 160 pieces

RT274-01-007

R3 775,45



Maico Eartip kit

- Suitable for the Maico EasyTymp, Maico touchTymp
- Contains eartips in 9 variations
- Suitable earpieces for adults and children
- Including equipment to remove the tips from the tympanometer
- Including ear loop to clean the eartips if needed
- 1 set contains 70 pieces

RT274-01-008

R3 775,45

- NB: specify



Audiology tools and consumables



RT274-02-009

HEINE SET OF 4 SPECULAE.
2.4mm, 3.0mm, 4.0mm, 5.0 mm
R994,75



RT274-02-001

OTOLIGHT
R517,50



RT274-02-004

REUSABLE OTOLIGHT TIPS X2
R256,45



RT274-07-048

PURESIL IMPRESSION
MATERIAL 1:1
R1159,20



RT274-07-049

PURESIL IMPRESSION MATERIAL
R1159,20



RT274-07-050

PURESIL IMPRESSION CATALYST
R454,25



RT274-07-043

PURETONE OTOSTOPS
SMALL
R359,95



RT274-07-044

PURETONE OTOSTOPS
MEDIUM
R359,95



RT274-07-045

PURETONE OTOSTOPS
LARGE
R359,95

Note: product supplied may not look exactly like the picture (EG: colour may differ depending on supplier stock)



Audiology tools and consumables



RT274-01-013

EEG ABRASIVE SKIN
PREPPING GEL
R465,75



RT274-01-012

EEG TEN 20 ELECTRODE PASTE
R2 202,25



RT274-01-012/2

SONAVELLE ELECTRODE GEL
R2 202,25



RT274-01-108

MAICO - TIP-TRODE, ADULT
GOLD FOIL 13MM (PKG OF 20)
R2 000,00



RT274-01-109

MAICO - TIP-TRODE, CHILD
GOLD FOIL 10MM (PKG OF 20)
R2 000,00



RT274-03-013

HARTMAN CROCODILE FORCEPS
R966,00



RT274-01-011/2

INSERT FOAM EARTIPS
SMALL-9.4MM (50 PER PACK)
SPECIFY SIZE WANTED
R859,05



RT274-01-011

INSERT FOAM EARTIPS
MEDIUM-13.7MM (50 PER PACK)
SPECIFY SIZE WANTED
R859,05



RT274-02-012

AUDIOMETER EARPHONE
COVERS, ACOUSTICALLY
TRANSPARENT. PACK OF 50.
R466,90

Note: product supplied may not look exactly like the picture (EG: colour may differ depending on supplier stock)



Audiology tools and consumables

RT274-01-081

PATH PROBE TIP "PT-A" (ADULT)
TRANSPARENT & PT-S
(NEWBORN) WHITE (BAG OF 2)

R465,80

RT274-01-020

1X REUSABLE GOLD
CUP ELECTRODE LEAD
ADULT-10MM

R448,50

RT274-01-021

1X REUSABLE GOLD
CUP ELECTRODE LEAD
PEADIATRIC-6MM

R448,50

RT274-01-030

PATH ELECTRODE SHIELDED
CABLE FOR ABR/ASSR

R8 590,00

RT274-01-082

PATH EP-DP OAE PROBE FOR
TEOAE/DPOAE WITH RED STRAIN RELIEF

R16 900,00

RT274-01-083

PATH ADVANCED POWER
SUPPLY / CHARGER ADA

R2 900,50

RT274-01-065

DISPOSABLE ELECTRODE
BLUE SENSOR (PACK OF 25)

R820,00

RT274-01-066

SNAP ELECTRODES
(BAG OF 10 STRIPS
X 3 ELECTRODES)

R650,00

RT274-01-067

SNAP ELECTRODES
(BAG OF 10 STRIPS
X 3 ELECTRODES)

R650,00

Note: product supplied may not look exactly like the picture (EG: colour may differ depending on supplier stock)



Audiology tools and consumables



RT274-01-084

PATH ADVANCED POWER
SUPPLY / CHARGER ADAPTER - EU
R2 900,50



RT274-01-001

PRE-CUT CLEANING FLOSS KIT
FOR OAE MACHINES &
TYMPANOMETER PROBES
R281,75



RT274-01-115

PATH EP-DP OAE PROBE FOR
TEOAE/DPOAE WITH RED
STRAIN RELIEF
R17 150,50



RT274-01-107

TIP TRODE ELECTRODE CABLE
FOR ECOCGH
R3 540,00



RT274-01-121

REPLACEMENT SET OF
HEARING AID SCREW DRIVERS
R402,00



RT274-02-002

BRUSH TOOL FOR REMOVING
WAX FROM HEARING AIDS
R143,75



RT274-01-119

REPLACEMENT
TUBE EXPANDER
R1 845,75



RT274-01-123

REPLACEMENT
TUBE THREADER
R258,75



RT274-01-131

CUT-OUT CUP, EAR SYRINGING,
PLASTIC, RECEIVER MAXIMUM 500ML
R517,50

Note: product supplied may not look exactly like the picture (EG: colour may differ depending on supplier stock)




Audiology tools and consumables



RT274-01-122

REPLACEMENT NEEDLE
NOSE PLIER
R460,00



RT274-01-143

CURETTES, PLASTIC
CERUMEN, PAEDIATRIC
R552,00



RT274-05-072

BRUSH TOOL FOR REMOVING
WAX FROM HEARING AIDS
R143,75



RT274-01-132

EAR IRRIGATION TIPS
R1782,50




RT274-01-120

REPLACEMENT SCALPEL
R322,00



RT274-07-046

SYRINGE, 10-15CM IN LENGTH
R1104,00



RT274-05-071

MAGNET TOOL FOR
REMOVING BATTERY
R161,00

Note: product supplied may not look exactly like the picture (EG: colour may differ depending on supplier stock)



Hearing Aid Care



RT274-07-021

PURETONE THIN-
WALLED TUBING
R399,05



RT274-07-022

PURETONE MEDIUM-
WALLED TUBING
R399,05



RT274-07-023

PURETONE THICK-
WALLED TUBING
R399,05



RT274-02-014

PURETONE BLOWER
R281,75



RT274-02-008

STETOCIP
R937,25



RT274-07-036

PAEDIATRIC HEARING AID CLIP
R327,75



RT274-07-013

BATTERY TESTER
R238,05



RT274-03-007

SYRINGE, EAR IRRIGATION, PLASTIC
R189,75



RT274-07-040

DRY AID, ELECTRONIC, FOR
HEARING AIDS, UNIVERSAL, KIT
R2 340,25

Note: product supplied may not look exactly like the picture (EG: colour may differ depending on supplier stock)

Audiology tool cleaning kit's



EARMOULD IMPRESSION STARTER KIT

CONTAINS: PURETONE UTILITY BAG, LED OTOLIGHT, PURESIL II SILICONE IMPRESSION MATERIAL (TUB WITH TUBE OF CATALYST), EAR MOULD SYRINGE, OTOBLOCKS (SMALL, MEDIUM, LARGE)

R4 496,50

RT274-07-055



PURETONE ADULT HEARING AID CARE KIT

CONTAINS: PURETONE UTILITY BAG, BATTERY TESTER, PURECLEAN CLEANING CUP, CLEANING TABLETS, STETOCLIP, PURETONE BLOWER, HEARING AID RETAINER, BTE HEARING AID RETAINER WITH CLIP (ADULT)

R1591.15

RT274-07-034



HEARING AID RETUBING KIT

CONTAINS: PURETONE UTILITY BAG, PURETONE BLOWER, NEEDLE NOSE PLIER, SCALPEL, SCREW DRIVER, TUBE EXPANDER, TUBE THREADER, CURVED TWEEZERS, SUPER GLUE

R3 680,00

RT274-01-118



Bone conduction hearing aids

Spec sheet:



PURETONE SOFTBAND MONAURAL BONE CONDUCTION HEARING AID

For mild-moderate conductive loss
Softband colour options: Black, blue, red or white
(Output 30-40)
6 Batteries (Size 13)
R15,410.00

**RT274-04-004
OFFER 2**

Spec sheet:



PURETONE DIGITAL HEADBAND MONAURAL BONE CONDUCTION HEARING AID

6 BATTERIES (SIZE 13)
R16 675

**RT274-04-005
OFFER 2**



Note: product supplied may not look exactly like the picture (EG: colour may differ depending on supplier stock)

Audiology
equipment and
consumables



Lebone Medical Supplies

Trading as

Amtronix Diagnostics

Partnering with you.



Front desk

Name. Thompho Lumadi

Phone. 012 054 5100/ 073 823 2380

Email. hear4u@amtronix.co.za

Audiology: Support, equipment advice, training

Name. Jade McEtywa

Phone. 081 046 1626

Email. audiochat@amtronix.co.za

Sales, quotes & orders

Name. Merle Pemberton

Phone. 078 350 5483

Email. diagnostic@amtronix.co.za

Technical support

Name. Tiaan Schoeman

Phone. 060 558 6199

Email. 911@amtronix.co.za

Arrangement of calibrations, call-outs & repairs

Name: Sorette Olivier

Phone. 012 054 5100/ 073 823 2380

Email. allears@amtronix.co.za



Website



Most **up-to date**
Catalogue of audiology
equipment on tender